

ATCAY Antibody



PACO07894

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ATCAY

Uniprot

Q86WG3

Synonyms:

ataxia, cerebellar, Cayman type; ATCAY; BNIP-H; CLAC; KIAA1872 ;

Immunogen:

Human ATCAY.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A