

# ATF4 Antibody



PACO07898

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ATF4

**Uniprot**

P18848

**Synonyms:**

activating transcription factor 4 (tax-responsive enhancer element B67); ATF4; CREB-2; CREB2; TAXREB67; TXREB ;

**Immunogen:**

Human ATF4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A