

# ATG12 Antibody



PACO07900

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

ATG12

**Uniprot**

O94817

**Synonyms:**

ATG12 autophagy related 12 homolog (S. cerevisiae); ATG12; APG12; APG12L; FBR93; HAPG12 ;

**Immunogen:**

Human ATG12.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A