

ATG16L1 Antibody



PACO07901

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ATG16L1

Uniprot

Q676U5

Synonyms:

ATG16 autophagy related 16-like 1 (*S. cerevisiae*); ATG16L1; APG16L; ATG16L; FLJ00045; FLJ10035; FLJ10828; FLJ22677; IBD10; WDR30 ;

Immunogen:

Human ATG16L1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A