

ATG5 Antibody



PACO07906

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ATG5

Uniprot

Q9H1Y0

Synonyms:

ATG5 autophagy related 5 homolog (S. cerevisiae); ATG5; APG5; APG5-LIKE; APG5L; ASP; hAPG5 ;

Immunogen:

Human ATG5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A