

ATP1B3 Antibody



PACO07918

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ATP1B3

Uniprot

P54709

Synonyms:

ATPase, Na⁺/K⁺ transporting, beta 3 polypeptide; ATP1B3; ATPB-3; CD298; FLJ29027 ;

Immunogen:

Human ATP1B3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A