

ATP2B2 Antibody



PACO07920

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ATP2B2

Uniprot

Q01814

Synonyms:

ATPase, Ca++ transporting, plasma membrane 2; ATP2B2; PMCA2; PMCA2a; PMCA2i ;

Immunogen:

Human ATP2B2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A