

ATP5A1 Antibody



PACO07923

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

ATP5A1

Uniprot

P25705

Synonyms:

ATP synthase, H⁺ transporting, mitochondrial F1 complex, alpha subunit 1, cardiac muscle; ATP5A1; ATP5A; ATP5AL2; ATPM; MOM2; OMR; ORM; hATP1 ;

Immunogen:

Human ATP5A1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A