

# ATP5L Antibody



PACO07929

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

ATP5L

**Uniprot**

O75964

**Synonyms:**

ATP synthase, H<sup>+</sup> transporting, mitochondrial F0 complex, subunit G; ATP5L; ATP5JG ;

**Immunogen:**

Human ATP5L.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A