

ATP5O Antibody



PACO07930

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ATP5O

Uniprot

P48047

Synonyms:

ATP synthase, H+ transporting, mitochondrial F1 complex, O subunit; ATP5O; ATPO; OSCP ;

Immunogen:

Human ATP5O.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A