

ATP6V1A Antibody



PACO07933

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ATP6V1A

Uniprot

P38606

Synonyms:

ATPase, H+ transporting, lysosomal 70kDa, V1 subunit A; ATP6V1A; ATP6A1; ATP6V1A1; HO68; VA68; VPP2; Vma1 ;

Immunogen:

Human ATP6V1A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A