

ATP6V1E1 Antibody



PACO07939

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ATP6V1E1

Uniprot

P36543

Synonyms:

ATPase, H+ transporting, lysosomal 31kDa, V1 subunit E1; ATP6V1E1; ATP6E; ATP6E2; ATP6V1E; P31; Vma4 ;

Immunogen:

Human ATP6V1E1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A