

AURKA Antibody



PACO07955

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

AURKA

Uniprot

O14965

Synonyms:

aurora kinase A; AURKA; AIK; ARK1; AURA; AURORA2; BTAK; MGC34538; STK15; STK6; STK7 ;

Immunogen:

Human AURKA.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A