

# AVIL Antibody



PACO07956

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

AVIL

**Uniprot**

O75366

**Synonyms:**

advillin; AVIL; ADVIL; DKFZp779O1812; DOC6; FLJ12386; MGC133244; p92 ;

**Immunogen:**

Human AVIL.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A