

AVPI1 Antibody



PACO07957

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

AVPI1

Uniprot

Q5T686

Synonyms:

arginine vasopressin-induced 1; AVPI1; PP5395; RP11-548K23.7; VIP32; VIT32 ;

Immunogen:

Human AVPI1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A