

AZIN1 Antibody



PACO07958

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

AZIN1

Uniprot

O14977

Synonyms:

antizyme inhibitor 1; AZIN1; MGC3832; MGC691; OAZI; OAZIN; ODC1L ;

Immunogen:

Human AZIN1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A