

BACE1 Antibody



PACO07970

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

BACE1

Uniprot

P56817

Synonyms:

beta-site APP-cleaving enzyme 1; BACE1; ASP2; BACE; FLJ90568; HSPC104; KIAA1149 ;

Immunogen:

Human BACE1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A