

BANP Antibody



PACO07978

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

BANP

Uniprot

Q8N9N5

Synonyms:

BTG3 associated nuclear protein; BANP; BEND1; DKFZp761H172; FLJ10177; FLJ20538; SMAR1; SMARBP1 ;

Immunogen:

Human BANP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A