

BCAP29 Antibody



PACO07983

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

BCAP29

Uniprot

Q9UHQ4

Synonyms:

B-cell receptor-associated protein 29; BCAP29; BAP29; DKFZp686M2086 ;

Immunogen:

Human BCAP29.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A