

BCAR1 Antibody



PACO07985

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

BCAR1

Uniprot

P56945

Synonyms:

breast cancer anti-estrogen resistance 1; BCAR1; CAS; CAS1; CASS1; CRKAS; P130Cas ;

Immunogen:

Human BCAR1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A