

BCAT1 Antibody



PACO07987

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

BCAT1

Uniprot

P54687

Synonyms:

branched chain aminotransferase 1, cytosolic; BCAT1; BCT1; DKFZp686E12175; ECA39; MECA39; PNAS-121; PP18 ;

Immunogen:

Human BCAT1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A