

BCL10 Antibody



PACO07990

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

BCL10

Uniprot

O95999

Synonyms:

B-cell CLL/lymphoma 10; BCL10; CARMEN; CIPER; CLAP; c-E10; mE10 ;

Immunogen:

Human BCL10.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A