

BCL2L10 Antibody



PACO07994

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

BCL2L10

Uniprot

Q9HD36

Synonyms:

BCL2-like 10 (apoptosis facilitator); BCL2L10; BCL-B; Boo; Diva; MGC129810; MGC129811 ;

Immunogen:

Human BCL2L10.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A