

BEX2 Antibody



PACO08006

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

BEX2

Uniprot

Q9BXY8

Synonyms:

brain expressed X-linked 2; BEX2; BEX1; DJ79P11.1 ;

Immunogen:

Human BEX2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A