

BHLHE41 Antibody



PACO08010

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

BHLHE41

Uniprot

Q9C0J9

Synonyms:

basic helix-loop-helix family, member e41; BHLHE41; BHLHB3; DEC2; SHARP-1; SHARP1 ;

Immunogen:

Human BHLHE41.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A