

BIN2 Antibody



PACO08014

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

BIN2

Uniprot

Q9UBW5

Synonyms:

bridging integrator 2; BIN2; BRAP-1 ;

Immunogen:

Human BIN2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A