

BLNK Antibody



PACO08022

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

BLNK

Uniprot

Q8WV28

Synonyms:

B-cell linker; BLNK; BASH; BLNK-S; LY57; MGC111051; SLP-65; SLP65 ;

Immunogen:

Human BLNK.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A