

# BMF Antibody



PACO08025

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

BMF

**Uniprot**

Q96LC9

**Synonyms:**

BCL2 modifying factor; Bmf; AW260063; AW411872 ;

**Immunogen:**

Human Bmf.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A