

# BMI1 Antibody



PACO08026

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

BMI1

**Uniprot**

P35226

**Synonyms:**

BMI1 polycomb ring finger oncogene; BMI1; MGC12685; PCGF4; RNF51 ;

**Immunogen:**

Human BMI1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A