

# BNIP1 Antibody



PACO08032

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

BNIP1

**Uniprot**

Q12981

**Synonyms:**

BCL2/adenovirus E1B 19kDa interacting protein 1; BNIP1; NIP1; SEC20; TRG-8 ;

**Immunogen:**

Human BNIP1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A