

BOLA1 Antibody



PACO08034

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

BOLA1

Uniprot

Q9Y3E2

Synonyms:

bolA homolog 1 (E. coli); BOLA1; CGI-143; MGC75015; RP11-196G18.18 ;

Immunogen:

Human BOLA1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A