

BTG4 Antibody



PACO08058

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

BTG4

Uniprot

Q9NY30

Synonyms:

B-cell translocation gene 4; BTG4; MGC33003; PC3B ;

Immunogen:

Human BTG4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A