

BTN3A3 Antibody



PACO08062

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

BTN3A3

Uniprot

O00478

Synonyms:

butyrophilin, subfamily 3, member A3; BTN3A3; BTF3 ;

Immunogen:

Human BTN3A3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A