

# MCMBP Antibody



PACO08070

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

MCMBP

**Uniprot**

Q9BTE3

**Synonyms:**

chromosome 10 open reading frame 119; C10orf119; FLJ13081; FLJ36756; MCM-BP; MCMBP ;

**Immunogen:**

Human C10orf119.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A