

ANAPC15 Antibody



PACO08072

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ANAPC15

Uniprot

P60006

Synonyms:

chromosome 11 open reading frame 51; C11orf51; DKFZp564M082; HSPC020 ;

Immunogen:

Human C11orf51.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A