

PACO08073

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

AAMDC

**Uniprot**

Q9H7C9

**Synonyms:**

chromosome 11 open reading frame 67; C11orf67; CK067; FLJ21035; MGC3367; PTD015 ;

**Immunogen:**

Human C11orf67.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A