

C16orf70 Antibody



PACO08082

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

C16orf70

Uniprot

Q9BSU1

Synonyms:

chromosome 16 open reading frame 70; C16orf70; C16orf6; FLJ12076; LIN10; lin-10 ;

Immunogen:

Human C16orf70.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A