

# C18orf8 Antibody



PACO08085

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RMC1

**Uniprot**

Q96DM3

**Synonyms:**

chromosome 18 open reading frame 8; C18orf8; FLJ23453; HsT2591; MIC1 ;

**Immunogen:**

Human C18orf8.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A