

C19orf43 Antibody



PACO08087

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

TRIR

Uniprot

Q9BQ61

Synonyms:

chromosome 19 open reading frame 43; C19orf43; MGC2803; fSAP18 ;

Immunogen:

Human C19orf43.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A