

C1orf106 Antibody



PACO08090

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

INAVA

Uniprot

Q3KP66

Synonyms:

chromosome 1 open reading frame 106; C1orf106; FLJ10901; MGC125608 ;

Immunogen:

Human C1orf106.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A