

# C1orf50 Antibody



PACO08092

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

C1orf50

**Uniprot**

Q9BV19

**Synonyms:**

chromosome 1 open reading frame 50; C1orf50; MGC955 ;

**Immunogen:**

Human C1orf50.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A