

PACO08094

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CHTOP

**Uniprot**

Q9Y3Y2

**Synonyms:**

chromosome 1 open reading frame 77; C1orf77; DKFZp547E1010; MGC131924; MGC86949; RP1-178F15.2; pp7704 ;

**Immunogen:**

Human C1orf77.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A