

C1S Antibody



PACO08100

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

C1S

Uniprot

P09871

Synonyms:

complement component 1, s subcomponent; C1S; FLJ44757 ;

Immunogen:

Human C1S.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A