

TMEM230 Antibody



PACO08102

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

TMEM230

Uniprot

Q96A57

Synonyms:

chromosome 20 open reading frame 30; C20orf30; HSPC274; dJ1116H23.2.1 ;

Immunogen:

Human C20orf30.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A