

C2orf49 Antibody



PACO08106

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

C2orf49

Uniprot

Q9BVC5

Synonyms:

chromosome 2 open reading frame 49; C2orf49; FLJ45759; MGC5509 ;

Immunogen:

Human C2orf49.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A