

# C3orf20 Antibody



PACO08107

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

C3orf20

**Uniprot**

Q8ND61

**Synonyms:**

chromosome 3 open reading frame 20; C3orf20; DKFZp434N1817 ;

**Immunogen:**

Human C3orf20.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A