

C4BPA Antibody



PACO08108

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

C4BPA

Uniprot

P04003

Synonyms:

597aa,67kd;

Immunogen:

Human C4BPA.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

Product Images

N/A

N/A