

# HAUS3 Antibody



PACO08110

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

HAUS3

**Uniprot**

Q68CZ6

**Synonyms:**

chromosome 4 open reading frame 15; C4orf15; DKFZp686l1868; IT1; MGC4701 ;

**Immunogen:**

Human C4orf15.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A