

# C5AR1 Antibody



PACO08111

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

C5AR1

**Uniprot**

P21730

**Synonyms:**

complement component 5a receptor 1; C5AR1; C5A; C5AR; C5R1; CD88 ;

**Immunogen:**

Human C5AR1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A