

TRAPPC13 Antibody



PACO08112

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

TRAPPC13

Uniprot

A5PLN9

Synonyms:

chromosome 5 open reading frame 44; C5orf44; FLJ13611; FLJ26957; MGC48585 ;

Immunogen:

Human C5orf44.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A